

## SpectrumViewer Frequently asked questions

Q What spectral data formats can SpectrumViewer read and write

A At the present time, SpectrumViewer can currently read and write the following spectral data formats

- o ORTEC data files: (\*.chn, \*.spc, and \*.spe, also known as IAEA ASCII)
- o Canberra CAM files: (\*.cnf, please be sure to see the next question for details)
- o IAEA ASCII data files: (\*.spe)
- o GADRAS data files: (\*.pcf and \*.asc)
- o IMS data files: (\*.phd)
- o IEC ASCII data files (\*.iec)
- o IEEE ASCII data files (\*.iee)
- o ANSI N42.42 XML data files: (\*.n42)
- o ASCII text files: (\*.txt, one or two column text preceded by a one line header in the following format: nChannels, liveTime, realTime, zero, gain, quad, sampleID string)

In addition, SpectrumViewer can also read the following two formats:

- o GR-135 binary data files (\*.dat)
- o AmpTek data files (\*.mca)

Q Is it possible to read and save spectral data collected in Canberra CAM (\*.cnf) format

A Yes, but you will need to copy the following dll files from your "GENIE2K.EXEFILES" (or equivalent) folder to the SpectrumViewer bin folder: Cam32.dll, ipc.dll, OS2WIN32.dll, Pcam.dll, and Sad.dll.

The files G2K\_VB.dll and SadVB.dll (along with several others) should already be there.